



LUSEAUX®

HYDRION® QT-10 TEST PAPER

MEASURES QUATERNARY AMMONIUM CONCENTRATIONS

GENERAL DESCRIPTION: A quick, reliable, easy-to-use, way to test quaternary ammonium-based sanitizers in food-service establishments, as required by federal, state, and local food safety codes. QT-10 is formulated for 2-chain quat compounds with measurements in the range of 0-400 parts per million (ppm). Each kit consists of a 15 ft. x 5/32 in. roll of premium test paper in a handy Hydrion® dispenser with a matching color chart. Approximately 100 tests per kit.

APPLICATION: Open dispenser and remove paper roll from the foil wrapping. Place the roll back into the dispenser and replace dispenser top. Tear a strip of HYDRION® QT TEST PAPER from the dispenser and dip into solution to be tested. Match at once to the color chart. Refer to sanitizer label for information on required concentrations.

Always consult your local health department for more information regarding sanitizing regulations for your area.

LUSEAUX® Laboratories, Inc.

DETERGENTS • CLEANERS • SANITIZERS

Gardena, CA 90247 • 800-266-1555

Kingman, AZ 86401 • 800-895-4583